The future is digital -

so are we!

PHYWE Systeme GmbH & Co. KG

Robert-Bosch-Breite 10 37079 Göttingen/Germany

P. +49 (0) 551 604 - 0 F. +49 (0) 551 604 - 107

info@phywe.com www.phywe.com



Scope of supply

1 x Cobra SMARTlink

- 1 x Charging cable type C with charging adapter
- 2 x Temperature probe type K
- 2 x Banana plug, stackable, red with 4 mm plug, 50 cm length
- 1 x Banana plug, stackable, blue with 4 mm plug, 50 cm length
- 1 x Manual (German, English)
- 1 x High quality packing carton incl. storage box

Technical data

- 10.1" LCD touch screen (16:10). WUXGA resolution 1200 x 1920
- Octa-core CPU 2.0 GHz, 4 GB RAM Main memory DDR3
- Internal memory 64 GB, expandable to 128 GB with TF card (max. 64 GB)
- Bluetooth 5.0, GPS, WiFi standard 802.11 a/b/g/n/ac
- Battery: Li-Po 10.000mAh
- Power supply: AC-DC adapter 100...240 V~/ 50...60 Hz. 5V / 2A output
- Weight: 1200 g (incl. packaging)

Confirmity

PHYWE Systeme GmbH & Co.KG hereby declares that the radio system type 12999-99 complies with the 2014/53/EU directive. The complete text of the EC Declaration of Conformity is available at: www.phvwe.com/en/ec-declaration

Safety instructions

/!



- Protect the device from dust, high humidity and vapours. Clean the device with a slightly damp, lint-free cloth. Strong cleaning agents or solvents are not suitable.
- Ensure that no liquid enters the openings in the housing.
- Do not open the device.
- The power input is protected by a self-resetting fuse. However, avoid short-circuit currents exceeding 20A. This can lead to the destruction of the device.

Disposal



- The packaging mainly consists of environmentally-friendly materials that should be returned to the local recycling stations.
- Do not dispose of this product with normal household waste.
- If this unit needs to be disposed of, please return it to the address that is stated on the back for proper disposal.

Manual

Cobra SMARTlink

Simple. Digital. Experimentation.



Cobra **PHYWE SMARTlink**

Thank you for choosing Cobra SMARTlink!

Cobra SMARTlink is a state-of-the-art, high-performance **3-in-1 tablet**. It is not only a great **measuring and display device** for digital science teaching, but also a full featured Android tablet.

Due to the brilliant 10.1 inch display combined with the powerful hardware and the **internal measuring sensors**, it is perfectly equipped for the demands of modern, digital science teaching.

Especially in student experiments **Cobra SMARTlink** is very good, easy and reliable to use. Besides the internal sensors, more than 40 further sensors (Cobra SMARTsense and Cobra DigiCart) via Bluetooth 🚯 and USB 🚓 connectable.



Internal sensors and much more! (**O**) $\langle \bullet \rangle$ (**"**11 2x 🌡 $\langle \mathbf{A} \rangle$ $\langle \mathbf{O} angle$ $\langle \widehat{\mathbf{n}} \rangle$ **Voltage** Speakers • Right/Left Temperature: thermocouple type K Acceleration Headphones · 3.5 mm socket • Measuring range: -30 V...+30 V Measuring range: -200...1200 °C Resolution: 20 mV Resolution: 0.4 °C Accuracy: +/- 1 % Accuracy: +/- 3 % (Microphone Data rate: 10 Hz • Data rate: 1 kHz (ight) Volume USB Typ C (A) Current Measuring range: 1-128 kLx $\langle \mathbf{Q} \rangle$ GPS Charging via USB Measuring range: -1A...+1A Resolution: 0.5 mA Resolution: 1 Lx Connection to PC Accuracy: +/- 1 % Accuracy: +/- 4 % Data rate: 10 Hz • Data rate: 1 kHz SD card slot $\langle \bullet \rangle$ 0n/Off Switch Ri (typical): 220 mΩ Humidity • Max. perm. Current: ± 1.1 A Memory expansion Measuring range: 0...100 % 🛞 uv Micro SD card Resolution: 0.1 % Measuring range: 0...400 W/m² Accuracy: +/- 4 % (0-80 %) Resolution: 0.1 W/m² Data rate: 10 Hz Accuracy: +/- 4 % • Data rate: 10 Hz



Camera • Front 2 MP Back 5 MP







Keyboard connection

Use

- Press the on/off switch () for more than 3 seconds to turn on the device
- Choose an app to record measurement data
- measureAPP for measurement with Cobra SMARTsense
- DigiCartAPP for measurement with Cobra DigiCart
- Press the on/off switch () for more than 3 seconds to turn off the device

Charging the battery

Before using the device for the first time the battery must be fully charged. Use the enclosed USB charging adapter and connect the Cobra SMARTlink with the USB-C cable.

Use of sensors

All internal sensors are automatically recognized by the pre-installed app measureAPP. Furthermore, there are over 40 Cobra SMARTsense sensors and Cobra DigiCart to connect to Cobra SMARTlink.

The sensors can be connected either wirelessly via Bluetooth (3) or directly via USB $\bullet \hookrightarrow$ to one of the 4 USB ports $\bullet \hookrightarrow$.