

# Viscosity of Newtonian and non-Newtonian liquids (rotary viscometer)



Physics

Mechanics

Mechanics of liquids &amp; gases

Chemistry

Physical chemistry

Viscosity

Applied Science

Engineering

Applied Mechanics

Fluid dynamics &  
Aerodynamics

Difficulty level

medium



Group size

2



Preparation time

45+ minutes



Execution time

45+ minutes


This content can also be found online at:



<http://localhost:1337/c/60046f73f3052e0003c4a4d6>

PHYWE

# General information



## Application

PHYWE



Fig 1:  
Experimental  
set-up

The difference in viscosity of Newtonian and Non-Newtonian fluids has huge implications on the behaviour of such fluids. As such the understanding of such differences is very important in fields such as fluid dynamics.

## Other information (1/2)

PHYWE

**Prior****knowledge****Main****principle**

The prior knowledge for this experiment is found in the Theory section.

The viscosity of liquids can be determined with a rotation viscometer, in which a motor with variable rotation speed drives a cylinder immersed in the liquid to be investigated with a spiral spring. The viscosity of the liquid generates a moment of rotation at the cylinder which can be measured with the aid of the torsion of the spiral spring and read on a scale.

## Other information (2/2)

PHYWE

**Learning  
objective****Tasks**

The goal of this experiment is to investigate the differences in behaviour between Newtonian and Non-Newtonian fluids.

1. Determine the gradient of the rotational velocity as a function of the torsional shearing stress for two Newtonian liquids (glycerine, liquid paraffin).
2. Investigate the temperature dependence of the viscosity of Castor oil and glycerine.
3. Determine the flow curve for a non-Newtonian liquid (chocolate).

## Theory (1/6)

PHYWE

If a liquid is between two plates and a force  $F$  acts along the plate in the direction of the  $x$  axis, the plate moves with velocity  $v$ . For Newtonian liquids the corresponding component of the shearing stress  $\tau$ .

$$\tau = \frac{F}{A} \quad (1)$$

is linked with the velocity gradient  $\frac{dv}{dx}$  as follows:

$$\tau = \eta \frac{dv}{dx} \quad (2)$$

( $\eta$  is the viscosity of the liquid and  $A$  the area of contact between the plate and the liquid.)

A number of substances (suspensions, emulsions) show a complex correlation between  $T$  and the integral velocity gradient  $D$  (non-Newtonian liquids). Hysteresis is also possible.

## Theory (2/6)

PHYWE




Fig. 2: Velocity gradient and shearing stress.




Fig. 3: Viscous and plastic flow of different substances 1. Newtonian (pure viscous) liquid 2. Dilatant liquid 3. Pseudoplastic liquid 4. Bingham (pure plastic) liquid 5./6. Quasiplastic liquid.

## Theory (3/6)

PHYWE

### Rotary viscometer

A rotary viscometer consists of an inner and an outer cylinder. The liquid to be investigated is located between them. At low rotational velocity the moment of rotation which is exerted on a cylindrical layer of liquid with a radius  $r$  and a height  $h$  conforms to the following relationship as a result of the rotation of the outer or inner cylinder:

$$T(r) = \tau \cdot 2\pi r h \cdot r \quad (3)$$

The shearing stress can be expressed by the measurable moment of rotation:

$$\tau(r) = \frac{T}{2\pi r^2 h} \quad (4)$$

In this case, the velocity gradient  $D$  is as follows:

$$D(r) = r \frac{d\omega}{dr} \quad (5) \quad \omega \text{ is the angular velocity.}$$

## Theory (4/6)

PHYWE

For Newtonian liquids eqn. (2) or eqn. (3) can be substituted in eqn. 1. Integration with the following limiting conditions:

$$\omega = 0 \text{ for } r = R_1$$

$$\omega = f \text{ for } r = R_2$$

( $R_1$  and  $R_2$  are the radii of the two cylinders) gives the following relationship between the measured moment of rotation and the angular velocity:

$$T = \frac{4\pi R_1^2 R_2^2 h}{R_2^2 - R_1^2} \eta f = C \eta f \quad (6)$$

where  $C$  is a device constant.

The above expression must be further corrected due to edge effects so that  $C$  becomes an empirical constant.

It is customary to use the average shearing stress acting on the surface of the two cylinders (2), which is obtained from the geometric or arithmetic mean of the shearing stresses ( $\tau_{mg}$  or  $\tau_{ma}$ )

$$\tau_{ma} = T \frac{R_1^2 + R_2^2}{4\pi h R_1^2 - R_2^2} \quad (7)$$

$$\text{or } \tau_{ma} = T \frac{1}{2\pi R_1 R_2 h} \quad (8)$$

## Theory (5/6)


PHYWE

Using Expression (4) the following is obtained for D

$$D_{\text{ma}} = \frac{R_2^2 + R_1^2}{R_2^2 - R_1^2} \cdot f \quad (9)$$

$$\text{or } D_{\text{mg}} = \frac{2R_1 + R_2}{R_2^2 - R_1^2} \cdot f \quad (10)$$

Fig. 4: Moment of rotation as a function of the frequency for a Newtonian liquid + Glycerine o Liquid paraffin



## Theory (6/6)

PHYWE


For non-Newtonian liquids T is no longer directly proportional to f nor is  $\tau$  proportional to D. There is an approximation formula which describes a relationship between T and  $\tau$  and between D and f.

For many liquids, the viscosity changes exponentially with the temperature  $T_{\text{abs}}$ :

$$\eta = Ae^{b/T_{\text{abs}}} \quad (\text{Andrade}) \quad (11)$$

$$\text{or } \log \eta = \frac{T_{\text{abs}} + b}{T_{\text{abs}} + c} \quad (\text{Vogel}) \quad (12)$$

Fig. 5: Moment of rotation as a function of frequency for a non-Newtonian liquid (chocolate at 302 K).




## Equipment

Position	Material	Item No.	Quantity
1	<a href="#">Magnetic stirrer with heater MRHei-Tec</a>	35752-93	1
2	<a href="#">Rotational viscometer</a>	18224-99	1
3	<a href="#">Right angle boss-head clamp</a>	37697-00	2
4	<a href="#">Supp.rod stainl.st.,50cm,M10-thr.</a>	02022-20	1
5	<a href="#">Spring balance holder</a>	03065-20	1
6	<a href="#">Support rod with hole, stainless steel, 10 cm</a>	02036-01	1
7	<a href="#">Magnetic stirring bar 30 mm, cylindrical</a>	46299-02	1
8	<a href="#">Separator for magnetic bars</a>	35680-03	1
9	<a href="#">Beaker, Borosilicate, low form, 600 ml</a>	46056-00	3
10	<a href="#">Beaker, Borosilicate, tall form, 250 ml</a>	46027-00	2
11	<a href="#">Glass rod, boro 3.3, l=200mm, d=5mm</a>	40485-03	1
12	<a href="#">Glycerol, 250 ml</a>	30084-25	2
13	<a href="#">Liquid paraffin, thick, 250 ml</a>	30180-25	1
14	<a href="#">Castor oil 250 ml</a>	31799-27	2
15	<a href="#">Acetone, chemical pure, 250 ml</a>	30004-25	3

PHYWE

# Setup and Procedure



## Setup

PHYWE

The experimental set-up is presented in Fig. 1. The rotary viscometer must be adjusted until it is exactly vertical. Use the adjustment screws to do this: they are located on the base of the support stand. There is a box level on the viscometer which allows one to check the exactness of the set-up's adjustment.

Screw the rotary cylinder on carefully (left-handed threads). Subsequently, lower the viscometer until the surface of the liquid exactly reaches the calibration mark on the rotary body in each case.



## Procedure (1/2)

PHYWE

Stir low viscosity fluids while heating to the desired measuring temperature with the aid of a magnetic stirrer and a magnetic stirring rod to rapidly achieve a uniform heat distribution. The temperature should always be measured in the immediate vicinity of the immersion cylinders. After the experimental temperature has been reached, turn off the heater. The temperature should remain constant for several minutes before measurements are begun, as the immersion cylinder must be in thermal equilibrium with the liquid. When thermal equilibrium has been reached, switch off the magnetic stirrer and determine the viscosity of the liquid. Since the moments of rotation which are measured in this experiment are very small, it is necessary to study the operating instructions of the rotary viscosimeter carefully and to follow them exactly.

After the measurement has been made, always clean the bar of the viscometer and the rotary cylinder carefully with water or acetone.

For glycerine and liquid paraffin, determine the dependence of the moment of rotation on the frequency in the range between 0.1 Hz and 1.0 Hz.

## Procedure (2/2)


PHYWE

For glycerine and castor oil, determine the dependence of the viscosity on the frequency in the temperature range between approximately 290 K and 350 K. For chocolate, determine the dependence of the moment of rotation of the frequency in the range between 0.1 Hz and 1.0 Hz at a temperature of approximately 303 K.

Other substances which are appropriate for investigation are: Newtonian substances: oils, ethylene glycol, etc.

Non-Newtonian liquids: paints (structural colour, hammer effect enamel), syrup, lubricants, chocolate spread, etc.

PHYWE




# Evaluation

## Results

PHYWE

At 303 K  
the  
viscosity of  
glycerine  
was  
calculated  
to be  
 $\eta = 680 \text{ cP}$



←  
Fig. 6:  
Temperature  
dependence of  
the viscosity of  
glycerine.

→  
Fig. 7:  
Temperature  
dependence of  
the viscosity of  
castor oil.

