





Velocity of sound using Kundt's tube and digital function generator



The goal of this experiment is to measure the speed of sound in air.


Physics	Acoustics	Sound generation & propagation	
Applied Science	Engineering	Applied Mechanics	Fluidynamics & Aerodynamics
 Difficulty level	 Group size	 Preparation time	 Execution time
medium	2	45+ minutes	45+ minutes

This content can also be found online at:



<http://localhost:1337/c/607d93a72934cc00036c4044>

PHYWE



General information

Application

PHYWE



Experimental setup

The speed of sound refers to the propagation speed of vibrations in a medium. As such the knowledge of this specific speed is very valuable in fields such as solid state physics and fluid dynamics.

This experiment investigates the specific speed of sound in air.

Other information (1/2)

PHYWE

Prior knowledge



No prior knowledge is required.

Main principle



Cork dust in a glass tube is set into tiniest motion by a sound wave. If the frequency of the sound wave matches the natural frequency of the volume in the glass tube, a standing wave will form. The cork dust then assembles in visible patterns that show the nodes of pressure and motion of the standing wave. From the length of the volume and the number of the nodes, the velocity of sound in the tube can be calculated for each natural frequency.

Other information (2/2)

PHYWE

Learning objective



The goal of this experiment is to measure the speed of sound in air.

Tasks



1. Determine the velocity of sound in air using Kundt's tube at different lengths of volume.

Theory (1/2)

PHYWE

The distance between two nodes of a standing wave is half its wavelength in case of parallel propagation of the waves, which can be assumed for this experiment. Therefore, the wave-length for each natural frequency can easily be determined and the frequency can be read directly from the digital function generator's display. Keeping in mind that the period T is the inverse of frequency

$$T = 1/\nu \quad (1)$$

and the wavelength is obtained from

$$\lambda = c/\nu = c \cdot T \quad (2)$$

The speed of sound can be found as the slope of the plot of wavelength vs. period.

Theory (2/2)

PHYWE

Alternatively, the wavelength can be calculated from the total number of nodes n and the length of the volume. For standing waves with one open end we obtain

$$\lambda = \frac{4L}{2n+1} \quad (3)$$

For two open ends one finds the following wavelength

$$\lambda = \frac{2L}{n} \quad (4)$$

Sound propagates in gases with varying speed depending on the gas temperature T_{gas} . From the equation of state for ideal gases one can derive the following relation

$$c \propto \sqrt{T_{gas}}$$


Comparison of sound speeds for different temperatures T_{gas} and T'_{gas} leads to

$$c = c' \cdot \sqrt{\frac{T_{gas}}{T'_{gas}}} \quad (5)$$

Equipment

Position	Material	Item No.	Quantity
1	PHYWE Digital Function Generator, USB	13654-99	1
2	Loudspeaker / Sound head, 8 ohms	03524-01	1
3	Tripod base PHYWE	02002-55	3
4	Universal clamp	37715-01	2
5	Kundt's apparatus	03475-88	1
6	Cork dust, 3 g	03477-00	1
7	Measuring tape, l = 2 m	09936-00	1
8	Thermometer -10...+50 °C	38034-00	1
9	Connecting cord, 32 A, 500 mm, red	07361-01	1
10	Connecting cord, 32 A, 500 mm, blue	07361-04	1
11	Rubber stopper,d=38/31mm,w/o hole	39260-00	1

PHYWE



Setup and Procedure

Setup

PHYWE

Perform the experimental set-up according to Fig. 1. Use the charging strip to distribute the cork dust over the glass tube. Use the piston to close one end of the glass tube. At the other end, set up the sound head in such manner, that the center of the sound head is in line with the center of the tube's cross sectional area. This way most of the sound will enter the tube and the amplification will stay bearable.

Connect the output of the digital function generator directly to the sound head via the connecting cords. Set the digital frequency generator to an amplitude of 3V to 5V and choose a sinusoidal signal and signal-type "out".



Fig. 1: Experimental setup

Procedure

PHYWE


Determine the velocity of sound in air. In order to do that, count the number of nodes for different natural frequencies for an open ended as well as a closed volume. For better statistics use the piston to examine at least two different lengths of volume.

To find the natural frequencies, start from 0 Hz and go up to maximum 4000 Hz in steps of 100 Hz.

When first patterns start to show, increase the frequency more slowly in steps of 10 Hz until the pattern does not wiggle anymore. When the whole pattern throughout the glass tube is absolutely static you found the natural frequency and generated a standing wave. Note down the number of nodes and the frequency before you continue to the next natural frequency. For each volume length find at least three natural frequencies for the open ended and the closed standing wave respectively.

Measure the room temperature during the experiment and note it down for later correction of the results.

PHYWE



Evaluation

Evaluation (1/2)

In the following, the evaluation of the obtained values is described with the help of example values. Your results may vary from those presented here.

Measurements with the plugged tube for two different volume lengths are shown in Table 1.

From Fig. 2, Fig. 3 and equation (2), we obtain the sound speed with 336 m/s and 330 m/s, respectively. The temperature during measurements was 19°C.

In the literature, a value of $c = 343$ m/s is given for room temperature of 20°C. Utilizing equation (5) the experimental velocity of sound at 20°C is obtained as 345 m/s and 339 m/s respectively.

L = 615 mm L = 350 mm

ν [Hz]	n	ν [Hz]	n
990	3	1150	2
1250	4	1600	3
2000	7	2250	4
2250	8	2700	5
2500	9	3200	6
3000	11	3654	7
3750	13		

Tab. 1: Measured values for different tube lengths L.

Evaluation (2/2)




Fig. 2: Measurements for L=615 mm.




Fig. 3: Measurements for L=350 mm.