

Charging curve of a capacitor / charging and discharging of a capacitor



Physics

Electricity & Magnetism

Electric current & its effects

Physics

Electricity & Magnetism

Simple circuits, resistors & capacitors

Applied Science

Engineering

Electrical Engineering

Properties of Electrical Devices

Applied Science

Engineering

Renewable Energy

Basic Principles



Difficulty level

easy



Group size

-



Preparation time

20 minutes



Execution time


30 minutes

This content can also be found online at:



<http://localhost:1337/c/64abba68c78b0a0002325fc1>

PHYWE



General information

Application

PHYWE



One example for the use of capacitors are radios

Capacitors are essential components in many everyday devices as they can store and discharge electrical energy. They are found in electronic circuits, motors, audio equipment, and power supply networks. A specific example is a computer, where capacitors are used in the power supply circuit to smooth the current and filter out disturbances. They are also present in radios, where they serve to tune into specific frequencies. Overall, capacitors contribute to the efficient and stable operation of our technology.

Other information (1/2)

PHYWE

**Prior****knowledge****Main****principle**

The prior knowledge required for this experiment is found in the theory section.

A capacitor is charged by way of a resistor. The current is measured as a function of time and the effects of capacitance, resistance and the voltage applied are determined.

Other information (2/2)

PHYWE

**Learning
objective****Tasks**

The goal of this experiment is to investigate the characteristics of capacitors.

To measure the charging current over time:

1. using different capacitance values C , with constant voltage U and constant resistance R .
2. using different resistance values (C and U constant)
3. using different voltages (R and C constant).


To determine the equation representing the current when a capacitor is being charged, from the values measured.

Theory (1/2)

PHYWE

The course of current with time, $I(t)$, when a capacitor C is charged through a resistor R at a fixed voltage U is determined from Kirchhoff's laws:

$$I(t) = \frac{U}{R} e^{-\frac{t}{RC}} \quad (1)$$




Capacitor charging circuit a) charging b) discharging

Theory (2/2)

PHYWE

If discharging curves are to be measured as well, the circuit as shown in the figure will be used.

Another experiment which could be carried out would be to determine unknown capacitance values from the charging and discharging curves with known resistance and charging function, or conversely to determine large resistance values at known capacitance.




Circuit for recording charging and discharging curves.

Equipment


Position	Material	Item No.	Quantity
1	Connection box	06000-00	2
2	Two-way switch, single pole	06005-00	1
3	Capacitor, 2 x 30 μ F	06007-00	1
4	Resistor 100 Ohm, 1W, G1	39104-63	1
5	Resistor 1 MOhm, 1W, G1	39104-52	4
6	bridge plug	06027-07	2
7	Capacitor 1 microF/ 100V, G2	39113-01	1
8	Capacitor 4,7microF/ 100V, G2	39113-03	1
9	PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1
10	Digital stopwatch, 24 h, 1/100 s and 1 s	24025-00	1
11	Digital multimeter, 600V AC/DC, 10A DC, 20M Ω	07124-12	1
12	Connecting cord, 32 A, 250 mm, red	07360-01	3
13	Connecting cord, 32 A, 250 mm, blue	07360-04	4

PHYWE


Setup and Procedure



Setup

PHYWE
excellence in science

- Set up the experiments as shown in the figures.
- Various resistance values R are established by series connection. The internal resistance of the digital multimeter and the setting time can be disregarded.
- R_1 is a protective resistor which limits the current when discharging (switch setting b).



Procedure

PHYWE

To measure the charging current over time:


1. using different capacitance values C , with constant voltage U and constant resistance R
2. using different resistance values (C and U constant)
3. using different voltages (R and C constant).

To determine the equation representing the current when a capacitor is being charged, from the values measured.



PHYWE

Evaluation



Evaluation (1/7)




Fig. 1: Course of current with time at different capacitance values; voltage and resistance are constant

$$U = 9V, R = 2,2M\Omega$$

The dependence of the current on the capacitance, the resistance and the voltage should be worked out from the measured values obtained by systematically varying the parameters.

1. First plot the measured values direct (Fig. 1) and then semilogarithmically (Fig. 2).




Fig 2: as Fig.1 but plotted semilogarithmically

Evaluation (2/7)




Fig. 1: Course of current with time at different capacitance values; voltage and resistance are constant

$$U = 9V, R = 2,2M\Omega$$

According to the figures the function takes the general form:

$$I(t) = I_0(U, R)e^{\alpha(U, R, C) \cdot t}$$

I_0 is not dependent on C as all curves begin at the same current values.





Fig 2: as Fig.1 but plotted semilogarithmically

Evaluation (3/7)

PHYWE

Fig. 3: α as a function of capacitance.

To investigate the dependence of the exponent on the capacitance, the slopes of the straight lines in Fig. 3 are plotted against capacitance, on a log-log basis.

A straight line with the slope $-0,98 \approx -1$ is obtained, so that:

$$I(t) = I_0(U, R)e^{\frac{\alpha'(U, R)}{C}t}$$

Evaluation (4/7)

PHYWE





Fig. 4: Course of current with time at different resistance values; capacitance and voltage are constant at $64\mu F$ and $9V$ respectively.

2. Straight lines with different slopes and different starting points are obtained. The dependence of the exponent on R is determined by plotting the log-log of the straight lines in Fig. 4.

A straight line with the slope -1.00 is obtained, so that

$$I(t) = I_0(U, R)e^{\frac{\alpha''(U)}{RC}t}.$$

Fig. 5: Exponent α as a function of resistance R .

Evaluation (5/7)

PHYWE




Fig. 6: Starting current I_0 of the measured values in Fig. 4 as a function of the resistance.

The straight line in Fig. 6 has a slope of $-0,99 \approx -1$ i.e.

$$I_0 = \frac{\beta(U)}{R}$$

Evaluation (6/7)

PHYWE




Fig. 7: course of current with time at various voltages
 $R = 4,4 M\Omega$, $C = 4 \mu F$.

3. All the straight lines in Fig. 7 have the same slope. The exponent is thus independent of the voltage U (this statement can also be made on the basis of dimensions).

The slope of the straight line is:

$$0,058 s^{-1} = \frac{1}{RC}$$

and therefore

$$RC = 17,24 s.$$

Evaluation (7/7)

PHYWE




Fig. 7: course of current with time at various voltages
 $R = 4,4M\Omega$, $C = 4\mu F$.

The starting current values I_0 for the measured values in Fig. 7 are plotted directly against the voltage values U in this case (Fig. 8).

A straight line with the slope $0,227 \frac{\mu A}{V} = \frac{1}{R}$

and therefore $R = 4,41M\Omega$

is obtained. Taken together, therefore, all the measured values give equation:

$$I(t) = \frac{U}{R} e^{(\frac{t}{RC})} \quad (1)$$