

## Magnetic moment in the magnetic field (Item No.: P2430400)



## Overview

## Short description

A conductor loop carrying a current in a uniform magnetic field experiences a torque. This is determined as a function of the radius, of the number of turns and the current in the conductor loop and of the strength of the external field.



#### **Student's Sheet**

Printed: 12.12.2017 11:23:46 | P2430400



Fig. 1: Experimental set-up for determining the torque due to a magnetic moment in the magnetic field.

### Equipment

Position No.	Material	Order No.	Quantity
1	Helmholtz coils, one pair	06960-00	1
2	Conductors, circular, set	06404-00	1
3	Torsion dynamometer, 0.01 N	02416-00	1
4	Coil holder for 02416-00	02416-02	1
5	Distributor	06024-00	1
6	PHYWE Power supply, universal, analogue display DC: 18 V, 5 A / AC: 15 V, 5 A	13503-93	1
7	PHYWE power supply, variable DC: 12 V, 5 A / AC: 15 V, 5 A	13540-93	1
8	Digital multimeter 2005	07129-00	2
9	Support base DEMO	02007-55	1
10	Support rod, stainless steel, 750 mm	02033-00	1
11	Right angle clamp expert	02054-00	2
12	Connecting cord, 32 A, 750 mm, red	07362-01	5
13	Connecting cord, 32 A, 750 mm, blue	07362-04	5

#### Tasks

Determination of the torque due to a magnetic moment in a uniform magnetic field, as a function

- 1. of the strength of the magnetic field,
- 2. of the angle between the magnetic field in the magnetic moment,
- 3. of the strength of the magnetic moment.



Robert-Bosch-Breite 10 D - 37079 Göttingen

**DHVWE** 

expert



## **Set-up and Procedure**

The experimental set-up is as shown in Fig. 1. Series connection is recommended so that the same magnetic field is induced in both coils. In the Helmholtz arrangement which can be built up with the spacing cross-members supplied, the coils are arranged reversed, so that the connections 1-1 or 2-2 should be joined (for series connection). In continuous operation the current in the Helmholtz coils should not exceed 3 A.

The connection wires to the coil carrier should hang loosely. They should be twisted together, so that no additional moment is produced.

The zero-point of the torsion balance should be checked frequently, since rapid rotary movements can displace the connecting leads.

Very small torques occur when measuring torque as a function of the Helmholtz coil current and of the angle. It is therefore recommended to use only the coil with 3 turns and to increase the coil current briefly (approx. To 6). The angles should be set at 15° intervals, by alternate use of the notches in the coil carrier.

## Theory and evaluation

With a closed conductor loop C , through which there flows a current I, a magnetic moment  $ec{m}$  is definded:

$$ec{m}=rac{I}{2}\oint_Cec{r} imes dec{r}=I\oint_Aec{d}\,\Omega$$
 .

A is any given area, the boundary of which is C . A magnetic field with flux density  $\dot{B}$  exerts a torque  $\dot{T}$  on a magnetic moment.

 $ec{T}=ec{m} imesec{B}$  (1)

If the magnetic field varies with position, the individual parts of the conductor loop are subjectet to different torques. It is therefore desirable to bring the conductor loop into a uniform magnetic field. Two coils, set up as shown in Fig. 1 and whose radius is equal to the distance between them (Helmholtz arrangement), are used to produce a uniform magnetic field.

For the present case, in which the conductor loop is a flat current ring with diameter d and n turns,

 $ec{m}=I\cdot n\cdotec{A}$  (2) $ec{m}ec{$ 

where  $\vec{A}$  is the vector of the area of the current ring. If a current I' flows in the Helmholtz coils, then, from (1):  $|\vec{T}| = c \cdot I \cdot n |\vec{A}| I' \cdot \sin \alpha$  (3)

where  $\alpha$  is the angle between  $\vec{B}$  and the plane vector  $\vec{A}$ , and c is a constant of these Helmholtz coils. The exponents of the different experimental set ups shown in table 1 are proving the above mentioned equations. Table 1



Robert-Bosch-Breite 10 D - 37079 Göttingen

#### **Student's Sheet**

Printed: 12.12.2017 11:23:46 | P2430400







From the regression line to the measured values of Fig. 2, with the exponential statement  $Y=A\cdot X^B$ 

the results are listed in Table 1.



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

info@phywe.de www.phywe.com

#### Student's Sheet Printed: 12.12.2017 11:23:46 | P2430400

TESS expert









Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

info@phywe.de www.phywe.com

# Student's Sheet

Printed: 12.12.2017 11:23:46 | P2430400









Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107 info@phywe.de www.phywe.com