
Ferromagnetism, paramagnetism and diamagnetism



The goal of this experiment is to investigate the different types of magnetism.

Physics	Electricity & Magne	etism		
Difficulty level	QQ Group size	D Preparation time	Execution time	
hard	2	45+ minutes	45+ minutes	







General information

Application







Magnetism is a class of physical phenomena that are mediated by magnetic fields. Electric currents and the magnetic moments of elementary particles give rise to a magnetic field, which acts on other currents and magnetic moments. Magnetism is one aspect of the combined phenomenon of electromagnetism.

There are three different types of magnetism: Ferromagnetism, paramagnetism and diamagnetism.







Robert-Bosch-Breite 10 37079 Göttingen

Theory



In a magnetic field, all of the materials are subject to more or less strong forces. Three types of materials can be distinguished.

Ferromagnetic materials: They are attracted by every magnet. Their permeability is very high, but not a constant. It depends on the external magnetic field intensity and the history of the material.

Paramagnetic materials: They are attracted into areas of high magnetic field intensity. Their relative permeability is slightly greater than one and it is independent of the field strength.

Diamagnetic materials: They are repelled from areas of high magnetic field intensity. Their relative permeability is slightly smaller than one and it is a material constant that is independent of the field strength.




Equipment

Position	Material	Item No.	Quantity
1	Tripod base PHYWE	02002-55	1
2	Support rod, stainless steel, 750 mm	02033-00	1
3	Right angle clamp expert	02054-00	1
4	Rod with hook	02051-00	1
5	U-magnet, large, U-shaped, limb length 130 mm, colored poles	06320-00	1
6	Pole pieces for iron core, U-shaped	06493-00	1
7	Sample set Ferro-Para-Diamagnetism, Ni, W, Bi	06340-00	1
8	Silk thread, I = 200 m	02412-00	1





Setup and Procedure



Procedure



- Fasten the thread with the nickel rod to the support system so that it is suspended above the space between the pole pieces.
- Observe its behaviour.
- Fasten the thread with the tungsten rod to the support system so that it is suspended between the pole pieces. Ensure that it can perform rotary oscillations.
- Observe the tungsten rod until it steadies.
- $\circ~$ Use the bismuth rod in the same manner as the tungsten rod.

Evaluation

Evaluation



The nickel rod is immediately attracted very strongly by the magnetic poles. It aligns itself in parallel with regard to the connecting line between the pole pieces.

In the end, the tungsten rod aligns itself in the direction of the connecting line between the pole pieces.

After it has stopped moving, the bismuth rod aligns itself perpendicularly with regard to the connecting line between the pole pieces.

Ferromagnetic and paramagnetic materials are attracted into the area of high magnetic field intensity, whereas diamagnetic materials are repelled.

