


# Hall effect in metals



Physics

Electricity & Magnetism

Magnetism & magnetic field

Applied Science

Engineering

Materials Science

Thermal & Electrical Properties

Applied Science

Engineering

Electrical Engineering

Properties of Electrical Devices



Difficulty level

hard



Group size

2



Preparation time

45+ minutes



Execution time

45+ minutes

This content can also be found online at:



<http://localhost:1337/c/600a96cc6080870003ab8b18>

PHYWE

# General information



## Application

PHYWE



Fig. 1: Experimental set-up for the Hall effect in metals.

The Hall effect is widely used everywhere where magnetic fields appear in solid matter. One especially important application is its use in semiconductor detectors.

This experiment offers the opportunity to gain a first understanding of the Hall effect in metals.

## Other information (1/2)

PHYWE

**Prior****knowledge****Main****principle**

The prior knowledge for this experiment is found in the Theory section.

The Hall effect in thin zinc and copper foils is studied and the Hall coefficient determined. The effect of temperature on the Hall voltage is investigated.

## Other information (2/2)

PHYWE

**Learning****objective****Tasks**

The goal of this experiment is to investigate the Hall effect in different metals.

1. The Hall voltage is measured in thin copper and zinc foils.
2. The Hall coefficient is determined from measurements of the current and the magnetic induction.
3. The temperature dependence of the Hall voltage is investigated on the copper sample.

## Theory (1/2)

PHYWE

If a current  $I$  flows through a strip conductor of thickness  $d$  and if the conductor is placed at right angles to a magnetic field  $B$ , the Lorentz force

$$\vec{F} = Q(\vec{v} \times \vec{B})$$

acts on the charge carriers in the conductor,  $n$  being the drift velocity of the charge carriers and  $Q$  the value of their charge. This leads to the charge carriers concentrating in the upper or lower regions of the conductor, according to their polarity, so that a voltage – the so-called Hall voltage  $U_H$  – is eventually set up between two points located one above the other in the strip:

$$U_H = \frac{R_H \cdot B \cdot I}{d}$$

$R_H$  is the Hall coefficient.

## Theory (2/2)

PHYWE

The type of charge carrier can be deduced from the sign of the Hall coefficient: a negative sign implies carriers with a negative charge (“normal Hall effect”), and a positive sign, carriers with a positive charge (“anomalous Hall effect”). In metals, both negative carriers, in the form of electrons, and positive carriers, in the form of defect electrons, can exist. The deciding factor for the occurrence of a Hall voltage is the difference in mobility of the charge carriers: a Hall voltage can arise only if the positive and negative charge carriers have different mobilities.

## Equipment

Position	Material	Item No.	Quantity
1	Hall effect, Cu, carrier board	11803-00	1
2	Hall effect, Zn, carrier board	11804-01	1
3	PHYWE Teslameter, digital	13610-93	1
4	Hall probe, tangential, protection cap	13610-02	1
5	PHYWE Universal measuring amplifier	13626-93	1
6	PHYWE Power supply, universal DC: 0...18 V, 0...5 A / AC: 2/4/6/8/10/12/15 V, 5 A	13504-93	1
7	PHYWE Variable transformer with digital display DC: 0...20 V, 12 A / AC: 0...25 V, 12 A	13542-93	1
8	Digital multimeter, 600V AC/DC, 10A AC/DC, 20 M $\Omega$ , 200 $\mu$ F, 20 kHz, -20°C... 760°C	07122-00	2
9	Coil, 300 turns	06513-01	2
10	Iron core, U-shaped, laminated	06501-00	1
11	Pair of pole pieces, plane, 30 x 30 x 48 mm	06489-00	1
12	Tripod base PHYWE	02002-55	1
13	Right angle clamp expert	02054-00	2
14	Universal clamp with joint	37716-00	1
15	Support rod, stainless steel, l = 250 mm, d = 10 mm	02031-00	1
16	On/off switch	06004-00	1
17	Connecting cord, 32 A, 750 mm, red	07362-01	7
18	Connecting cord, 32 A, 750 mm, blue	07362-04	5
19	Connecting cord, 32 A, 750 mm, black	07362-05	2

PHYWE



# Setup and Procedure

## Setup (1/2)

PHYWE

The layout follows Fig. 1 and the wiring diagram in Fig. 2.

- Arrange the field of measurement on the plate midway between the pole pieces.
- Carefully place Hall probe in the centre of the magnetic field.
- The measuring amplifier takes about 15 min. to settle down free from drift and should therefore be switched on correspondingly earlier.
- To keep interfering fields at a minimal level, make the connecting cords to the amplifier input as short as possible.
- Take the transverse current  $I$  for the Hall probe from the powersupply unit 13536.93. It can be up to 15 A for short periods.

## Setup (2/2)

PHYWE




Fig. 2 Circuit diagram for the Hall effect.

## Procedure (1/2)

PHYWE

The Hall probe will show a voltage at the Hall contacts even in the absence of a magnetic field, because these contacts are never exactly one above the other but only within manufacturing tolerances. Before measurements are made, this voltage must be compensated with the aid of the potentiometer as follows:

- Disconnect the transverse current  $I$ .
- Set the measuring amplifier to an output voltage of 1 V, for example, by adjusting the compensation-voltage. ( $R_e = 104 \Omega$ , amplification = 105).
- Connect the transverse current.
- Twist the connecting cords between hall voltage sockets and amplifier input in order to avoid as much as possible stray voltages..
- Adjust the compensating potentiometer, using a screwdriver, until the instrument again shows an output voltage of 1 V.

## Procedure (2/2)

PHYWE

- Repeat this operation several times to obtain a precise adjustment.

The determination of the Hall voltage is not quite simple since voltages in the microvolt range are concerned where the Hall voltages are superposed by parasitic voltages such as thermal voltages, induction voltages due to stray fields, etc. The following procedure is recommended:

- Set the transverse current  $I_x$  to the desired value.
- Set the field strength  $B$  to the desired value (on the power supply, universal, 13500.93).
- Set the output voltage of the measuring amplifier to about 1.5 V by adjusting the compensation-voltage.
- Using the mains switch on the power supply unit, switch the magnetic field on and off and read the Hall voltages at each on and off position of the switch (after the measuring amplifier and the multi-range meter have recovered from their peak values). The difference between the two values of the voltage, divided by the gain factor 105, is the Hall voltage  $U_H$  to be determined.

PHYWE



## Evaluation

## Results (1/5)

PHYWE

The measurements for copper shown in Fig. 3 are related by the expression  $U_H \sim B$ . Linear regression using the relation  $U_H = a + bB$  shows these values to be represented by a straight line with the slope  $b = -0.0384 \cdot 10^{-6} \text{ m}^2/\text{s}$  and a standard deviation  $s_b = 0.0004 \cdot 10^{-6} \text{ m}^2/\text{s}$ . From this, with  $d = 18 \cdot 10^{-6} \text{ m}$  and  $I = 10 \text{ A}$ , we derive the Hall coefficient

$$R_H = -(0.576 \pm 0.006) \cdot 10^{-10} \text{ m}^3/\text{As}$$




Fig. 3: Hall voltage as a function of magnetic induction  $B$ , using a copper sample.

## Results (2/5)

PHYWE

The measurements shown in Fig. 4 are related by the expression  $U_H \sim I$ . Linear regression using the relation  $U_H = a + bI$  shows these values to be represented by a straight line with the slope  $b = -0.770 \cdot 10^{-6} \text{ V/A}$  and a standard deviation  $s_b = 0.0005 \cdot 10^{-4} \text{ V/A}$ . From this, with  $d = 18 \cdot 10^{-6} \text{ m}$  and  $B = 0.25 \text{ T}$ , we derive the Hall coefficient

$$R_H = -(4.31 \pm 0.004) \cdot 10^{-10} \text{ m}^3/\text{As}$$




Fig. 4: Hall voltage as a function of current  $I$ , using a copper sample.

## Results (3/5)

PHYWE

The zinc sample shows an anomalous Hall effect, in that the Hall voltage has a positive sign. The measurements shown in Fig. 5 are related by the expression  $U_H \sim B$ . Linear regression using the relation  $U_H = a + bB$  shows these values to be represented by a straight line with the slope  $b = -20.7 \cdot 10^{-6} \text{ m}^2/\text{s}$  and a standard deviation  $s_b = 0.4 \cdot 10^{-6} \text{ m}^2/\text{s}$ . From this, with  $d = 25 \cdot 10^{-6} \text{ m}$  and  $I = 12 \text{ A}$ , we derive the Hall coefficient

$$R_H = -(4.31 \pm 0.01) \cdot 10^{-11} \text{ m}^3/\text{As}$$




Fig. 5: Hall voltage as a function of magnetic induction  $B$ , using a zinc sample.

## Results (4/5)

PHYWE

The measurements shown in Fig. 6 for zinc are related by the expression  $U_H \sim I$ . Linear regression using the relation  $U_H = a + bI$  shows these values to be represented by a straight line with the slope  $b = 0.40 \cdot 10^{-6} \text{ V/A}$  and a standard deviation  $s_b = 0.01 \cdot 10^{-4} \text{ V/A}$ . From this, with  $d = 25 \cdot 10^{-6} \text{ m}$  and  $B = 0.25 \text{ T}$ , we derive the Hall coefficient

$$R_H = -(4.00 \pm 0.01) \cdot 10^{-11} \text{ m}^3/\text{As}$$




Fig. 6: Hall voltage as a function of current  $I$ , using a zinc sample.

## Results (5/5)

PHYWE

If the sample temperature is varied, we find, disregarding disturbing thermal voltages, that the Hall voltage in metals is not temperature dependent.

If the measured values of the Hall coefficients are compared with those given in the literature

(copper:  $R_H = -0.53 \cdot 10^{-10} \text{ m}^3/\text{As}$ ; zinc:  $R_H = 10 \cdot 10^{-10} \text{ m}^3/\text{As}$ ), it is noteworthy that the values for zinc show a distinct difference. This, it may be presumed, is attributable to disturbing secondary effects particularly at the contact points. Among these we may mention the Ettinghausen effect, the Peltier effect, the Seebeck effect, the first Righi-Leduc effect and the first Ettinghausen-Nernst effect. It is possibly due also to impurities in the test material (99.95 % purity).